



- Packet of 2 x 15ml bottles of NJ Phillips Lubricant recommended for use with NJ Phillips instruments.
- Use after cleaning to lubricate seal rings, cylinders and valves, or when changing to new seal rings, cylinder or valves.



**njphillips.com**  
*Our expertise is in your hands.*

LOCKED BAG 8, GOSFORD NSW, AUSTRALIA 2250  
P +61 2 4340 2044 F +61 2 4340 1991 | E [njp1@njphillips.com.au](mailto:njp1@njphillips.com.au)