

# F-Grip Dial-a-Dose Injector Range

- For use with a wide variety of injectable solutions.
- Lightweight, accurate and easy to use.
- Superior comfort grip and easy-off needle nut.
- Supplied with 2 needles.
- Sterilisable by boiling.
- Feed tube included.



## F-GRIP DIAL-A-DOSE

MODEL	DOSE RANGE	INCREMENTS	CODE
1ml	0.1ml – 1ml	0.1ml	PAS1394
2ml	0.1ml – 2ml	0.2ml	PAS1396
5ml	1ml – 5ml	0.5ml	PAS1398

10ml available on request, minimum order applies.

## F-GRIP DIAL-A-DOSE BOTTLE ATTACHMENT

MODEL	DOSE RANGE	INCREMENTS	CODE
1ml	0.1ml – 1ml	0.1ml	PAS1395
2ml	0.1ml – 2ml	0.2ml	PAS1397
5ml	1ml – 5ml	0.5ml	PAS1399

10ml available on request, minimum order applies.

KIT NO.	PART NAME	COMPONENTS
WX1754	1ml Major Kit	Cylinder, valves, spring & seal ring.
WX1755	2ml Major Kit	Cylinder, valves, spring & seal ring.
WX1756	5ml Major Kit	Cylinder, valves, spring & seal ring.
WX1768	125ml Bottle Adaptor Assembly	
WX1809	250ml Bottle Adaptor Assembly	

PLEASE ORDER BY KIT AND PART NAME

## CONFIGURATIONS TO SUIT STANDARD LUER NUT



## CONFIGURATIONS TO SUIT MIM NUT



## ACCESSORIES

SH1613 SPANNER



Customise this instrument with your brand colour and logo. Contact us for details.



**njphillips.com**  
Our expertise is in your hands.

LOCKED BAG 8, GOSFORD NSW, AUSTRALIA 2250

P +61 2 4340 2044 F +61 2 4340 1991 | E njp1@njphillips.com.au